



Rev. Sept. 2022



Food Reactivity Data

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result	
Cereals	Rice powder	-	Seeds	Peanut	-	Fruits	Orange	-	
	Polished rice	-		Roasted almond	-		Kiwi	-	
	Sprouted brown rice	-		Roasted cashew nut	-		Apple	-	
	Glutinous rice	-		Roasted macadamia nut	-		Peach	-	
	Buckwheat	-		Roasted pistachio nut	-		Banana	-	
	Japanese barnyard millet	-		Roasted walnut	-	Coconut milk	-		
	Foxtail millet	-		Roasted pecan nut	-	Coffee bean	-		
	Common millet	-		White sesame	-	Cocoa powder	-		
	Wheat	-		Black sesame	-	Green tea	-		
	Rye	-	Sea foods	Salmon	-	Coffee* Cocoa*Tea	Roasted laver	-	
	Barley	-		Mackerel	-		Seaweeds	Hijiki (boiled)	-
	Malt	-		Squid	-			Brown seaweed (raw)	- *
	Oat	-		Octopus	-			Kelp 200-fold extract	-
	Adlay	-		Black tiger prawn	-	Spices		Cumin	-
	Corn	-		Snow crab	-		Poppy seeds	-	
Corn flour	-	Scallop adductor		-	Black pepper		-		
Amaranth	-	Pollack roe		-	White pepper		-		
Quinoa	-	Salmon roe		-	Red pepper		-		
Tubers* Starch	Sago palm powder	-		Meats	Beef		-	Ginger	-
	Yam	-	Chicken		-		Cinnamon	-	
	Potato starch	-	Pork		-		Coriander	-	
Beans	Soybean	-	Eggs		Chicken egg		+	Polysaccharide Thickener	Turmeric
	Red bean (azuki)	-		Processed chicken egg	+		Low methoxyl pectin 50-fold extract		-
	Common bean	+		Dried chicken egg	+	Guar gum 1000-fold extract	-		
	Concord pole	-	Dairy products	Milk	-	Others	Gelatin	-	
	Broad bean	-		Skim milk	-				
Green stuff	Fungus	Tomato		-					
		Onion		-					
		Spinach	-						
Mushroom	-								

Test Result	
-	Negative
+	Positive (reaction positive)

*The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test Pro II (using Recommended Extraction Method), except Kelp 200-fold extract, Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.