



Rev. Sept. 2022



Food Reactivity Data

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result
Cereals	Rice powder	-	Seeds	Peanut	+ *1	Fruits	Orange	-
	Polished rice	-		Roasted almond	-		Kiwi	-
	Sprouted brown rice	-		Roasted cashew nut	-		Apple	-
	Glutinous rice	-		Roasted macadamia nut	-		Peach	-
	Buckwheat	-		Roasted pistachio nut	-		Banana	-
	Japanese barnyard millet	-		Roasted walnut	-		Coconut milk	-
	Foxtail millet	-		Roasted pecan nut	-	Coffee bean	-	
	Common millet	-		White sesame	-	Cocoa powder	-	
	Wheat	-		Black sesame	-	Green tea	-	
	Rye	-		Salmon	-	Roasted laver	-	
	Barley	-	Mackerel	-	Hijiki (boiled)	-		
	Malt	-	Squid	-	Brown seaweed (raw)	-		
	Oat	-	Octopus	-	Kelp 200-fold extract	-		
	Adlay	- *2	Black tiger prawn	-	Cumin	-		
	Corn	-	Snow crab	-	Poppy seeds	-		
Corn flour	-	Scallop adductor	-	Black pepper	-			
Amaranth	-	Pollack roe	-	White pepper	-			
Quinoa	-	Salmon roe	-	Red pepper	-			
Tubers • Starch	Sago palm powder	-	Meats	Beef	-	Ginger	-	
	Yam	-		Chicken	-	Cinnamon	-	
	Potato starch	-		Pork	-	Coriander	-	
Beans	Soybean	-	Eggs	Chicken egg	-	Polysaccharide Thickener	Turmeric	-
	Red bean (azuki)	-		Processed chicken egg	-		Low methoxyl pectin 50-fold extract	-
	Common bean	-	Dried chicken egg	-	Guar gum 1000-fold extract	-		
	Concord pole	-	Dairy products	Milk	-	Others	Gelatin	-
	Broad bean	-		Skim milk	-			
Green stuff	Tomato	-	Fungus	Mushroom	-			
	Onion	-						
	Spinach	-						

Test Result	
-	Negative
+	Positive (reaction positive)

*1 The test result may be negative if it is included in high concentrations.

*2 The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test Pro II (using Recommended Extraction Method), except Kelp 200-fold extract, Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.