



## Food Reactivity Data

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result
Cereals	Rice powder	-	Sea foods	Salmon	-	Coffee • Cocoa • Tea	Coffee bean	-
	Polished rice	-		Mackerel	-		Cocoa powder	-
	Sprouted brown rice	-		Squid	-		Green tea	-
	Glutinous rice	-		Octopus	-	Seaweeds	Roasted laver	-
	Buckwheat	-		Black tiger prawn	+ *1		Brown seaweed	-
	Japanese barnyard millet	-		Snow crab	+		Cumin	-
	Foxtail millet	-		Scallop adductor	-		Poppy seeds	-
	Common millet	-		Pollack roe	-	Black pepper	-	
	Wheat	-		Salmon roe	-	White pepper	-	
	Rye	-		Beef	-	Red pepper	-	
	Barley	-	Chicken	-	Wasabi (green horseradish)	-		
	Malt	-	Pork	-	Ginger	-		
	Oat	-	Chicken egg	-	Curry powder (spice mix)	-		
	Adlay	- *2	Dried chicken egg	-	Cinnamon	-		
	Corn	-	Dairy products	Milk	-	Coriander	-	
	Corn flour	-		Skim milk	-	Turmeric	-	
Amaranth	-	Green stuff	Tomato	-	Polysaccharide Thickener	Low methoxyl pectin 50-fold extract	-	
Quinoa	-		Onion	-		Guar gum 1000-fold extract	-	
Tubers • Starch	Sago palm powder		-	Spinach	-	Others	Gelatin	-
	Yam	-	Fungus	Mushroom	-			
	Potato starch	-		Orange	-			
Beans	Soybean	-		Kiwi	-			
	Red bean (azuki)	-		Apple	-			
	Common bean	-		Peach	-			
	Concord pole	-		Banana	-			
	Broad bean	-		Coconut milk	-			
Seeds	Peanut	-						
	Roasted almond	-						
	Roasted cashew nut	-						
	Roasted macadamia nut	-						
	Roasted pistachio nut	-						
	Roasted walnut	-						
	White sesame	-						
Black sesame	-							

Test Result	
-	Negative
+	Positive (reaction positive)

N.T. : Not Tested

\*1 The test result may be negative if it is included in high concentrations.

\*2 The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test Pro II for Crustacean (using Recommended Extraction Method), except Kelp 200-fold extract, Polysaccharide and Thickener samples.