



Rev. Sept. 2022



Food Reactivity Data

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result
Cereals	Rice powder	-	Seeds	Peanut	-	Fruits	Orange	-
	Polished rice	-		Roasted almond	-		Kiwi	-
	Sprouted brown rice	-		Roasted cashew nut	-		Apple	-
	Glutinous rice	-		Roasted macadamia nut	-		Peach	-
	Buckwheat	-		Roasted pistachio nut	-		Banana	-
	Japanese barnyard millet	-		Roasted walnut	+		Coconut milk	-
	Foxtail millet	-		Roasted pecan nut	+	Coffee bean	-	
	Common millet	-		White sesame	-	Cocoa powder	-	
	Wheat	-	Black sesame	-	Green tea	-		
	Rye	-	Salmon	-	Roasted laver	-		
	Barley	-	Mackerel	-	Hijiki (boiled)	-		
	Malt	-	Squid	-	Brown seaweed (raw)	-		
	Oat	-	Octopus	-	Kelp 200-fold extract	-		
	Adlay	-*	Black tiger prawn	-	Spices	Cumin	-	
	Corn	-	Snow crab	-		Poppy seeds	-	
	Corn flour	-	Scallop adductor	-		Black pepper	-	
Amaranth	-	Pollack roe	-	White pepper		-		
Quinoa	-	Salmon roe	-	Red pepper		-		
Tubers • Starch	Sago palm powder	-	Meats	Beef		-	Ginger	-
	Yam	-		Chicken		-	Cinnamon	-
	Potato starch	-		Pork		-	Coriander	-
Beans	Soybean	-		Eggs		Chicken egg	-	Turmeric
	Red bean (azuki)	-	Processed chicken egg			-	Polysaccharide Thickener	Low methoxyl pectin 50-fold extract
	Common bean	-	Dried chicken egg		-	Guar gum 1000-fold extract		-
	Concord pole	-	Dairy products	Milk	-	Others	Gelatin	-
	Broad bean	-		Skim milk	-			
Green stuff	Fungus	Tomato		-				
		Onion		-				
		Spinach	-					
Mushroom	-							

Test Result	
-	Negative
+	Positive (reaction positive)

*The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test Pro II (using Recommended Extraction Method), except Kelp 200-fold extract, Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.